

BEFORE APPLYING ORTHOSIS:

- Check leg for any neurological or vascular changes.
- Check device to be sure that lock and lock safety are engaged.



TO APPLY ORTHOSIS:

- Unfasten or loosen all straps.
- Slide device onto foot and align tension springs with ankle joint (maleoli).
- Align struts with long bones of the calf, making sure tension springs remain aligned with ankle joint.
- Refasten and tighten straps in this order:
 - Top ankle strap / Calf strap / Foot cradle straps
- Straps should be snug enough to allow one finger to be inserted between strap and skin.



FOR DYNAMIC TENSION:

- Swing lock lever safety away from lock lever.
- Disengage lock lever by slightly extending foot.



TO INSPECT THE FIT:

- Check alignment and adjust straps for comfort.
- Check for any neurological or vascular changes.



TO REMOVE ORTHOSIS:

- Lock device into position by engaging lock lever safeties with lock levers.
- Unfasten/loosen straps and carefully slide orthosis off foot.





ADVANCE[®] DYNAMIC ROM - ANKLE ORTHOSIS

RESTORE MOTION, RETURN TO FUNCTION



Empi, Inc. | A DJO Global Company

T 800.328.2536 **D** 651.415.9000

205 Highway 22 East | Clear Lake, SD 57226 | U.S.A.

Empi.com

Together in Motion™

