

# LOWER EXTREMITY DEVICES

## Lower Extremity Measurement Chart – GL Devices

### Patient Information:

Name: \_\_\_\_\_

- \_\_\_\_\_ A      Circumference: 1 1/2" (3.8 cm) below Groin
- \_\_\_\_\_ B      Circumference: Largest point of Calf
- \_\_\_\_\_ C      Length: Groin to Knee Center (*not center of Knee Cap*)
- \_\_\_\_\_ D      Length: Knee Center (*not center of Knee Cap*) to Medial Malleolus
- \_\_\_\_\_ E      Length: Foot Length

Date of Birth: \_\_\_\_\_

Height: \_\_\_\_\_

Weight: \_\_\_\_\_

Device: \_\_\_\_\_

Side: Left      Right

Flexion      Extension

Range of Motion requirements to use the GL Knee Devices:

> GL Knee Flexion - Leg must be able to bend to at least 48 degrees of Flexion

> GL Knee Extension - Leg must be able to bend to at least 50 degrees of Flexion

